


40480 QA

BCTuning Check 2023-06-26 09:31:41

[Log-File](#)

Measurement Step #1


Reflection





Reflection

	Frequency (at Minimum)				Factor (at Minimum)				Phase (at System Freq.)			
	Value	Low Spec	High Spec	Unit	Value	Low Spec	High Spec	Unit	Value	Low Spec	High Spec	Unit
0 deg	123.04	122.80	123.30	MHz	0.18	0.00	0.30		27.69			deg
90 deg	123.14	122.80	123.30	MHz	0.14	0.00	0.30		32.75			deg
Delta	100.0	-140.0	100.0	kHz								

Transmission





Transmission (Maximum)

	Value	Low Spec	High Spec	Unit
Frequency	123.57	122.80	123.60	MHz
Decoupling	-24.29	-40.00	-16.00	dB

BCCS modes

Switching of different modes in BCCS successful.

Results are in specification

Program finished successfully