

40480 TuneUp

Coil Power Losses 2023-06-26 09:19:32

[Log-File](#)

Phantom Center

	Value	Low Spec	High Spec	Unit
Phantom Center X	-1.06	-20.00	20.00	mm
Phantom Center Z	3.17	-20.00	20.00	mm

Phantom Diameter

	Value	Low Spec	High Spec	Unit
Phantom Diameter X	169.53	150.00	190.00	mm
Phantom Diameter Z	167.61	150.00	190.00	mm

Name	Value	Low_Spec	High_Spec	Unit
Reference Power from TALES	1604.63	500.00	2600.00	W

Coil Power Losses result 3883.21 W saved

Coil power losses saved successfully!

Results are in specification

Program finished successfully